

→
REF. M2-MTD
MEDICINE TROLLEY



TECHNICAL FEATURES



Width:	67 cm
Depth:	50 cm
Height:	90 cm
Ø of the wheels:	10 cm
Weight:	42 kg
Max. Load:	50 kg / drawer
Material:	Aluminum alloy & ABS
Surface:	White ABS
Color:	Grey

SPECIFICATIONS



- Soft closing sliding rails
- Lockable wheels
- Central lock with key

INCLUDES



- 4 drawers
- 2 waste bins
- Needle box
- 304 Stainless steel wire basket
- Extendable side tray
- Dividers to organize the drawer

