

→  
REF. EL12  
ELECTRIC TILT TABLE



## TECHNICAL FEATURES



Length:	189 cm
Width:	64 cm
Height:	48-95 cm
Weight:	85 kg
Max. Load:	175 kg
Motors:	2
Foam thickness:	5 cm
Density:	36 kg/m <sup>3</sup>
Upholstery:	MED PVC ●

## MOTOR FUNCTIONS



- Height (flat position): 48-95 cm
- 0° to 85° inclination angle

## MANUAL FUNCTIONS



- Footrest adjustable
- Armrest adjustable

## INCLUDES



- Hand controller
- Medical control box
- Breather hole
- 4 wheels with brake system

## CERTIFICATIONS




- 93/42/EEC Medical Device Directive

# Cran

→  
TECHNICAL DATA SHEET

CONFY CARE

→  
REF. EL12  
ELECTRIC TILTTABLE



**KG**  
max 225 kg

