

→
REF. FIX-MT2
FIXED TABLE



TECHNICAL FEATURES



2-section table

Length:	195 cm
Width:	70 cm
Height:	75 cm
Weight:	28 kg
Max. Load:	225 kg
Foam thickness:	5 cm
Density:	36 kg/m ³
Upholstery:	MED PVC ●

MANUAL FUNCTIONS



- Backrest 75° inclination

INCLUDES



- Breather hole
- Paper roll holder
- Folding legs

CERTIFICATIONS



- 93/42/EEC Medical Device Directive

→
REF. FIX-MT2
FIXED TABLE



KG
max 225 kg

