



TECHNICAL DATA SHEET

REF. CFALO5F ALUMINIUM PORTABLE MASSAGE TABLE



TECHNICAL FEATURES

 \rightarrow

Product size (Length x width x height)

Open:	185x70x63-84 cm
Folded:	93x71x22 cm
Weight:	17 kg
Max. Load:	250 kg
Foam thickness:	6 cm
Upholstery:	PU ●

MANUAL FUNCTIONS

 \rightarrow

- Height (flat position): 63-84 cm
- Backrest 70° inclination

SPECIFICATIONS

 \Rightarrow

- Aluminium welded frame
- Round corners

INCLUDES



- Breather hole
- Hole headrest
- Arm support
- Carry bag

CERTIFICATIONS CE



- 93/42/EEC Medical Device Directive

