

NEW REF. 110004*

REF. EL032

ELECTRIC TREATMENT TABLE
CAMILLA DE TRATAMIENTO ELÉCTRICA
TABLE DE TRAITEMENT ÉLECTRIQUE
LETTINO ELETTRICO PER TRATTAMENTI
ELEKTRISCHE BEHANDLUNGSLIEGE



195 x 66 cm



49,5 - 93 cm



69 kg



Max.

250 kg



1



5 cm - 36 kg/m³



1100041911



1100042411



1100043421



1100043521



MOTOR FUNCTIONS/ FUNCIONES MOTORIZADAS/ FONCTIONS MOTORISÉES/ FUNZIONI MOTORIZZATE/ MOTORISIERTE FUNKTIONEN

- Height/ Altura/ Hauteur/ Altezza/ Höhe: 49,5-93 cm


MANUAL FUNCTIONS/ FUNCIONES MANUALES/ FONCTIONS MANUELLES/ FUNZIONI MANUALI/ MANUELLE FUNKTIONEN

- Adjustable armrests - Headrest (-45° to +40°) and Footrest (0° to +75°) inclination by gas lift
- Reposabrazos ajustables - Reposacabezas (-45° a +40°) y Reposapiés (0° a +75°) por pistón de gas
- Accoudoirs réglables - Repose-tête (-45° et +40°) et Repose-pied (0° et +75°) par vérin à gaz
- Braccioli regolabili - Poggiatesta (-45° a +40°) e Poggiapiedi (0° a +75°) pistone a gas
- Einstellbare Armlehne - Kopfstütze (-45° bis +40°) und Fußstütze (0° bis +75°) mittels Gasdruckfeder

INCLUDES/ INCLUYE/ INCLUS/ INCLUDE/ ENTHÄLT

- Peripheral switch for height
- Hole headrest
- 4 retractable big wheels system for easier movement
- Paper roll holder
- Barra perimetral para regular la altura
- Reposacabezas con agujero facial
- Sistema de 4 ruedas grandes retráctiles para facilitar el movimiento
- Portarrollos
- Barre périphérique pour régler la hauteur
- Repose-tête avec trou visage
- Système de 4 grandes roues rétractables pour faciliter sa mobilité
- Porte-rouleau
- Barra perimetrale per regolare l'altezza
- Poggiatesta con foro viso

NEW REF. 110004*
REF. EL032



- Sistema di 4 grandi ruote retrattili per un facile spostamento
- Portarotolo
- Peripherie-Schalter für Höhenverstellung
- Kopfstütze mit Gesichtslöch
- 4 großen, einziehbaren Räder um die Beweglichkeit zu erleichtern
- Papierrollenhalter

CERTIFICATIONS/ CERTIFICACIONES/ CERTIFICATIONS/ CERTIFICAZIONI/ ZERTIFIZIERUNGEN

(EU) 2017/745

- Class I Medical Device. European regulation of medical devices
- Dispositivo Médico Clase I. Reglamento Europeo de Productos Sanitarios
- Dispositif Médical de Classe I. Règlement Européen des Dispositifs Médicaux
- Dispositivo Medico di Classe I. Regolamento Europeo dei Dispositivi Medici
- Medizinprodukt der Klasse I. Europäische Verordnung über Medizinprodukte

