



FACE SHIELD

DISPOSABLE FACE SHIELD

WITH SOFT ELASTICATED STRAP

LIGHTWEIGHT FACE SHIELD

Lightweight face shield with an anti-fog coating protects the face and eyes from splashing or spraying of infectious or hazardous substances. The disposable face shields are manufactured from optically clear polyester for maximum visibility, whilst the soft elastic strap and PU foam headband provide additional, long-term comfort and durability.

The product can be worn with or without additional face protection, such as prescription or protective eyewear and face masks.

Our registered design face shield is manufactured using a patent-protected process.



KEY FEATURES

- ✦ Anti-fog protection
- ✦ Optically clear visor for maximum visibility
- ✦ PU foam headband for long-term comfort
- ✦ Soft elasticated strap
- ✦ Lightweight design
- ✦ Protects the face and eyes from splashing or spraying
- ✦ Can be worn with or without eyewear protection and face masks
- ✦ CE Category III

SPECIFICATIONS

SPEC	DETAIL
PRODUCT COLOUR	OPTICALLY CLEAR (NO TINT)
PRODUCT TYPE	FACE SHIELD
SIZE	330MM X 220MM
OPTICAL CLASS	1 (+/-0.06)
MECHANICAL STRENGTH	S
FIELD OF USE	DROPLETS OR SPLASH OF LIQUIDS
APPLICATION	MEDICAL, GENERAL
STANDARDS	BS EN166:2002 CERT REF: 2777/14610-01/E00-00 CATEGORY III – COMPLEX DESIGN
SHELF LIFE	2 YEARS
STORAGE CONDITIONS	STORE IN ORIGINAL PACKAGING OUT OF DIRECT SUNLIGHT AT 23°C AND 50% RH

MATERIAL DATA

MATERIAL	DESCRIPTION	SIZE
VISOR FILM	OPTICAL GRADE PET	0.175MM THICK
FOAM HEADBAND	SOFT GRADE PU	25MM THICK
FOAM HEADBAND ADHESIVE	PSA	NA
ELASTIC HEADBAND	MEDICAL GRADE WOVEN ELASTIC	19MM THICK

INSTRUCTIONS FOR USE

Carefully remove the protection film before use. To put on and take off products ensure the elastic is wrapped securely around the head with foam against the forehead. Only wear products of a suitable size that are correctly adjusted. Products which are loose may not stay in place and if too tight will be uncomfortable. The size of this product is available in one standard size with 340mm elastic headband.

ADDITIONAL SPECIFICATIONS

Medical grade materials / Single use, cleaning not recommended / Latex Free / Can be used for prolonged periods as it does not create distortion or fatigue / Visor suitable for use with prescription lenses and surgical masks / Always check before using / Manufactured in the UK

TECHNICAL DRAWINGS

